created  
to protect  
your life

# S1 MASK

- Ⓘ Istruzioni d'uso e manutenzione  
\_\_\_\_\_
- Ⓙ Use and maintenance instructions  
\_\_\_\_\_
- Ⓚ Instrucciones de uso y mantenimiento

**DG srl**

## S1 MASK

Mascherina chirurgica multiuso in tessuto antibatterico TYPE IIR

DIASTAR GROUP SRL  
Corso Allamano 34  
10095 Grugliasco Torino

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## S1 MASK

Multipurpose surgical mask in TYPE IIR antibacterial fabric

**DIASTAR GROUP SRL**  
Corso Allamano 34  
10095 Grugliasco - ITALY

**Class I medical device,  
non-sterile, without  
measuring function**

|                   |                              |
|-------------------|------------------------------|
| <b>S1 MASK</b>    | IIR type surgical mask       |
| <b>SIZE</b>       | Standard                     |
| <b>BE CAREFUL</b> | Wear covering mouth and nose |
| <b>REUSABLE</b>   | Up to 50 washings            |

Use and maintenance instructions

### 1. PRODUCT DESCRIPTION / RESISTANCE OF THE FABRIC TO THE SUBSTANCES

Personal protective clothing in antibacterial fabric

This instruction manual refers to clothing made with material that has particular characteristics and that gives the products made with it the following characteristics:

- **Water repellency**
- **Waterproofness**
- **Antibacteriability**
- **Breathability**
- **Non-permeability of contaminants such as viruses and bacteria**

These Class I Medical Devices allow for adequate protection for the user since they are designed to protect the entire working environment and also the wearer from splashes of potentially contaminated liquids.

The use of these garments must always be accompanied with the knowledge of the instructions for use and does not exclude the application of the required hygiene regulations. The garments do not have an expiration date, but their use causes an inevitable consumption due to mechanical agents, which is why the user must assess whether the integrity of the garment is sufficient to guarantee its protective efficiency. The garments are made in such a way as not to obstruct users' work movements and guarantee them the maintenance of an internal environment with a correct degree of humidity, allowing perspiration. They have also been tested and verified for surgical use (Linting and Particulate test), therefore they ensure the level of hygiene required by the UNI EN 13795-2: 2009 standard on the generation of particles in a dry state, the Linting and Particulates are well below the limit imposed by the standard. Therefore they can be used in operating rooms and in environments where this characteristic is required (laboratories of genetic microbiology, electronics etc ...).

### 2. USE OF THESE INSTRUCTIONS

Read very carefully

The above symbol indicates the need to increase the degree of attention, in particular when reading these indications which provide information and warnings for the correct use of the gown and for its maintenance. The instructions always accompany the product, even in a reduced form and with reference to the safety criteria when using it. The instructions must be respected and

enforced to guarantee safety at work in relation to the use of the masks. Failure to comply with the instructions or partial and inadequate compliance can expose the user and the environment to risks; it may also cause early deterioration of special repellency treatments.

### 3. WARRANTY CONDITIONS

The products we sell are covered by a guarantee regarding the structure and use of suitable materials, under the following conditions:

1) The warranty is valid for a period of six (6) months, is given on the seams, on the good manufacture of the fabric and the product.

**2) The guarantee on the barrier capacity to bacterial penetration is guaranteed for 5 washes.**

The washes are carried out in water at 40 ° C or at 60 ° C or in the presence of washing with hypochlorite for fabrics.

The first contact with highly aggressive acids such as caustic soda over 90% and hydrogen peroxide at 130% terminates the warranty period as the article in these particular cases protects the operator only once (except direct verification for reuse). Protection is given by the waterproof membrane.

The warranty does not include damage deriving from excessive stress such as for example the use of the product after an anomaly has been detected, e.g. tearing or excessive wear, the use of inadequate operating methods as well as failure to observe the instructions of use. Furthermore, it does not include improper use and the manufacturer assumes no responsibility for contamination of the operator, who must always be well informed on the use of this device and on its correct use in conjunction with other protective devices.

**NOTICE: If the use of the guarantee is deemed necessary, please indicate the following data:**

- 1) **Typology**
- 2) **Purchase date**
- 3) **Detailed description of the problem**
- 4) **Production batch number**

### 4. TECHNICAL ASSISTANCE

The manufacturer as indicated on the label remains available to the customer for any type of request.

### 5. CERTIFICATIONS:

CERTIFIED PRODUCT Dir 93/42 / CEE

### 6. SAFETY WARNINGS

ATTENTION: IF THE DEVICE IS DAMAGED, IT MUST NOT BE USED

CAREFULLY CHECK THAT ALL SEAMS ARE ALIGNED. STRETCH THE FABRIC TO THE END TO ENSURE THAT THERE ARE NO LACERATIONS OR HOLES. PULL THE ELASTICS BY APPLYING A SLIGHTLY HIGHER STRENGTH THAN WHICH IS OBTAINED BY WEARING THE GARMENT.

#### 6.1 General warnings

Read the instructions carefully before using the product.

Failure to observe the safety warnings can cause serious damage to the user and those around him.

Incorrect use will impair the performance of the product and void any warranty claims.

#### 6.2 Normal use

The internal part of the fabric is coated with silver particles, these carry out the biocidal action as demonstrated by the test and carried out according to the UNI EN ISO 20743: 2007 standard on Staphylococcus aureus (ATCC6538).

• **Wash preferably below 40 ° C to keep water repellency longer and to avoid cooking microorganisms on the fabric itself and stain it indelibly. Washing with normal detergent which is already activated at 30 ° C normally present on the market, on this particular type of fabric, is able to eliminate most of the stains.**

For washing, normal disinfectant additives for washing machines referring to fabrics can be used, which are normally found on the market (Napisan, whitening products, half a glass of bleach, etc.) as long as they are used in cold water. Do not use pure bleach to thoroughly wet the fabric, it will quickly lose its water repellency. Remember that the fabric itself has biocidal efficacy thanks to the silver particles fixed in it permanently.

• **If required by operational uses, garments can be autoclaved at 134 ° C.**

• The garment used in the healthcare sector must be subjected to the sterilization protocol envisaged by the Istituto Superiore della Sanità, therefore immediately after use in the healthcare environment, if it has been contaminated by organic substances or body fluids, it must be subjected to decontamination as required from

the DM September 28, 1990. The garment must therefore be immediately sprayed with disinfectants whose chemical composition has at least 2/3 of ethyl alcohol and 1/3 of isopropyl with quaternary ammoniums and activators of recognized efficacy HIV-HBV-HCV bactericidal, fungicidal, tuberculocidal; let the disinfectant act for 1 minute and wait 3-5 minutes for it to be perfectly dry, then pack it in a preferably sealed bag and take it to washing (if an external laundry is used, inform the same to handle with safety guards for the material potentially contaminated and wash separately from other garments). The clean garment, if required by surgical needs, must then be wrapped and sterilized in an autoclave suitable for guaranteeing the distribution of any spores present.

## 7. STORAGE AND MAINTENANCE

The material before use must be stored in the following ways

- Inside its casing
- Protected against external agents of any type which could alter its performance characteristics
- **Far from heat or flame sources**

## 8. TRANSPORT

The product must be transported in its original packaging or with containers that do not alter its characteristics, during transport it must be protected from dust or other elements that can alter its functional conditions.

## DATA SHEET

### Summary of UNI EN 14126: 2004 results

Resistance to penetration of contaminated liquids under hydrostatic pressure (ISO 16604). After 5 washing cycles at 40 ° C.

| Class                | 1                                   | 2        | 3       | 4     | 5      | 6      |
|----------------------|-------------------------------------|----------|---------|-------|--------|--------|
| Hydrostatic pressure | 0 kPa                               | 1,75 kPa | 3,5 kPa | 7 kPa | 14 kPa | 20 kPa |
| Results              | <b>20 kPa – Class 6 compliant</b> ▲ |          |         |       |        |        |

Resistance to the penetration of infectious agents due to mechanical contact with substances containing contaminated liquids (ISO 22610). After 5 washing cycles at 40 ° C

| Class          | 1                                     | 2       | 3       | 4       | 5       | 6    |
|----------------|---------------------------------------|---------|---------|---------|---------|------|
| Passage time t | ≤10min                                | 15<t≤30 | 30<t≤45 | 45<t≤60 | 60<t≤75 | t>75 |
| Results        | <b>t &gt;75 – Class 6 compliant</b> ▲ |         |         |         |         |      |

Resistance to the penetration of contaminated liquid aerosols (ISO 22611). After 5 washing cycles at 40 ° C.

| Class                 | 1                                       | 2           | 3       |
|-----------------------|---|-------------|---------|
| Penetration ratio log | 1 < log ≤ 3                             | 3 < log ≤ 5 | log > 5 |
| Results               | <b>log &gt; 5 – Class 3 compliant</b> ▲ |             |         |

Resistance to the penetration of contaminated solid particles (ISO 22612). After 5 washing cycles at 40 ° C.

| Class               | 1  | 2               | 3   |
|---------------------|--|-----------------|-----|
| Penetration log ufc | 2 < log ufc ≤ 3                          | 1 < log ufc ≤ 2 | ≤ 1 |
| Results             | <b>Log ufc 5 1 - Class 3 compliant</b> ▲ |                 |     |

SUMMARY OF RESULTS UNI EN 14325: 2005 - UNI EN 13034:2009  
Determination of dimensional variations (30 ° C) (UNI EN 25077: 1996)

| Requirements | max ± 3% fabric   max ± 5% mesh | Compliant |
|--------------|---------------------------------|-----------|
| Results      | length: -1.5 % / width: -0.5 %  |           |

Burst resistance (Uni EN ISO 13938-1). After 5 washing cycles at 40 ° C.

| Class                  | 1                                    | 2        | 3         | 4         | 5         | 6         |
|------------------------|--------------------------------------|----------|-----------|-----------|-----------|-----------|
| Perforation resistance | > 40 kPa                             | > 80 kPa | > 160 kPa | > 320 kPa | > 640 kPa | > 850 kPa |
| Results                | <b>445 kPa - Class 4 compliant</b> ▲ |          |           |           |           |           |

Abrasion resistance (UNI EN 530: 1996 method 2). After 5 washing cycles at 40 ° C.

| Class            | 1                                      | 2     | 3     | 4      | 5      | 6      |
|------------------|--|-------|-------|--------|--------|--------|
| Number of cycles | > 10                                   | > 100 | > 500 | > 1000 | > 1500 | > 2000 |
| Results          | <b>&gt; 2000 – Class 6 compliant</b> ▲ |       |       |        |        |        |

Tensile strength of the seams (ISO 13935-2: 2001). After 5 washing cycles at 40°C.

| Class         | 1                                  | 2      | 3      | 4       | 5       | 6       |
|---------------|------------------------------------|--------|--------|---------|---------|---------|
| Pulling force | > 30 N                             | > 50 N | > 75 N | > 125 N | > 300 N | > 500 N |
| Results       | <b>N 400 – Class 5 compliant</b> ▲ |        |        |         |         |         |

Resistance to penetration of liquids UNI EN 6530/05 (Classification en 14325). For penetration resistance class 2 must be achieved for at least one of the specified chemical reagents. For liquid repellency class 3 must be achieved for at least one of the chemical reagents. After 5 washing cycles at 40 ° C.

| Class                              | 1                                 | 2    | 3    |
|------------------------------------|-----------------------------------|------|------|
| Penetration                        | < 10%                             | < 5% | < 1% |
| H <sub>2</sub> SO <sub>4</sub> 30% | <b>0.3% - Class 3 compliant</b> ▲ |      |      |
| Na OH 10%                          | <b>0.0% - Class 3 compliant</b> ▲ |      |      |
| o-Xylene                           | 11,3 % - NC                       |      |      |
| Butan-1-ol                         | <b>5.2% - Class 1 compliant</b> ▲ |      |      |

| Class                              | 1                                  | 2     | 3     |
|------------------------------------|------------------------------------|-------|-------|
| Repellency                         | > 80%                              | > 90% | > 95% |
| H <sub>2</sub> SO <sub>4</sub> 30% | <b>99.0% - Class 3 compliant</b> ▲ |       |       |
| Na OH 10%                          | <b>99.6% - Class 3 compliant</b> ▲ |       |       |
| o-Xylene                           | 38,4 % - NC                        |       |       |
| Butan-1-ol                         | 58,1 % - NC                        |       |       |

Further tests:

R.d.P. 1000979 of the Centro Tessile Cottoniero S.p.A.: surgical drapes, coveralls and gowns for operating blocks. Dry particle generation (Linting). Product standard UNI EN 13795-2: 2009 test standard ISO 9073-10: 2003 + UNI EN 13795-2: 2009 Linting and Particulate, Coincent of Linting = 2.6 (4 max value); Index parti culate matter (1PM) = 2.4 (3.5 max value). USL Trento prof. N. 2009.0036802 test in operating environment (repeated analysis with R.d.P. 5801915 - confirms the excellent microbial reduction).

| Parameter   | Cotton medical gown |      | Gown made of non-absorbent fabric |      |
|---|---------------------|------|-----------------------------------|------|
|   | Clean               | Used | Clean                             | Used |
| Total microbial load (u.f.c./24 cm <sup>2</sup> )   | 6                   | 80   | 0,5                               | 4,5  |
| Staphylococcus aureus (u.f.c./24 cm <sup>2</sup> )  | 0                   | 2    | 0                                 | 1    |
| Pseudomonas aeruginosa (u.f.c./24 cm <sup>2</sup> ) | 0                   | 0    | 0                                 | 0    |
| Ifomiceti (u.f.c./24 cm <sup>2</sup> )              | 0,5                 | 1    | 0,5                               | 1,5  |

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